

UniFusion SP50 Syringe Pump

Color touch LCD display

Multiple infusion modes

Long battery time

High water proof level



4.3 " color touch LCD display with text and image

±2% high accuracy to make sure safe infusion

3 infusion modes to satisfy basic infusion requirements

Use 5, 10, 20, 30, 50/60ml syringes

Easy calibration and automatic recognition of syringe size

Programmable and support change rate during infusion

Waterproof level IP34

Up to 11 hours battery life

Inter-lockable and free combination between syringe pump and infusion pump