

UniFusion SP50 Vet Syringe Pump

Color touch LCD display

Multiple infusion modes

Long battery time

High water proof level



UniFusion SP50 Vet Syringe Pump works both on battery and AC/or DC power supply, which is perfectly adapted for veterinary use. Our multi-function syringe pump bring precision, easy of use, safety and versatility for your day to day fluid therapy requirements.

Uses any standard syringes sizes 5, 10, 20, 30, 50/60ml and stores up to 20 editable brands

4.3 inch color touch screen, parameter setting and editing directly

Multi-function operation with 3 modes:

Rate Mode, Time Mode, Body-weight mode

Micro Mode - friendly design for small pets

Flow rates range from 0.1ml/h to 1,500ml/h with increment of 0.01ml/h

Double CPU for increased safety, visual, audible, symbolized full range alarm system, IrDA, WIFI, Nurse Call, RS232, data export etc

History record more than 5000 entries

Standard: 5.5 hours; Optional: 11 hours (@25ml/h)

IP34 waterproof, handle included and inter-lockable design