

MedRena

UniFusion VP50 Pro Infusion Pump

Color touch LCD display

Multiple infusion modes

Long battery time

High water proof level



4.3 " color touch LCD display with text and image

±5% high accuracy to make sure safe infusion

8 infusion modes to satisfy different infusion requirements

Double pressure sensor and double air bubble sensor reassure infusion safety

Programmable and support change rate during infusion

Big history records and drug library storage

Open system and close system optional

Automatic and manual bolus

DERS (Drug error reduction system)

DPS (Dynamic pressure system)

Waterproof level IP34

Up to 9 hours battery life

Inter-lockable and free combination between syringe pump and infusion pump