

UniFusion VP50 Vet Infusion Pump

Color touch LCD display

Multiple infusion modes

Long battery time

High water proof level



UniFusion VP50 Vet Infusion Pump works both on battery and AC/or DC power supply, which is perfectly adapted for veterinary use. Our multi-function infusion pump bring precision, easy of use, safety and versatility for your day to day fluid therapy requirements.

Uses any standard IV sets and stores up to 20 editable brands, close system available

4.3 inch color touch screen, parameter setting and editing directly

Multi-function operation with 4 modes:

Rate Mode, Time Mode, Body-weight mode and Drip mode

Micro Mode - friendly design for small pets

Flow rates range from 0.1ml/h to 1,500ml/h with increment of 0.01ml/h

Double CPU for increased safety, Ultrasonic air-in-line detector

IrDA, WIFI, Nurse Call, RS232, data export etc

Historic record more than 5000 entries

Standard: 4.5 hours; Optional: 9 hours (@25ml/h)

IP34 waterproof, handle included and inter-lockable design